

Automation of Sample Preparation in Ion Chromatography

Praha, 2nd October 2014

Ing. Roman Repáš
ACO Sales Manager (IC and SP)

The world leader in serving science

Agenda

- **Thermo Scientific Dionex Sample Preparation Products**
- **Automation of Sample Preparation in Ion Chromatography**

pH and Conductivity Measurement


Inline Filtration

Inline Sample Preparation

Matrix Elimination using Dionex InGuard Chemistries

Matrix Elimination with Selective Analyte Trapping


Thermo Scientific Dionex Sample Prep Product Line



Novel & Innovative Solutions

ThermoFisher
SCIENTIFIC

How Does Accelerated Solvent Extraction Work?



AutoTrace 280 : Automated SPE Extraction of large volume of liquid samples

- Automated SPE instrument for large volume aqueous samples (20 mL to 4 L)
- Extracts 6 samples with unattended operation
- All SPE parameters automatically controlled
- EPA approved for drinking water extraction



ThermoFisher
SCIENTIFIC

Dionex SolEx Solid-Phase Extraction Cartridges

- Cartridges designed for fast, and easy preconcentration of contaminants prior to analysis
- Offered in a variety of chemistries and sizes to fit your sample extraction needs
 - Silica-Based SPE Cartridges
 - SolEx C8
 - SolEx C18
 - SolEx C8-Clean (Phthalate-Free)
 - SolEx C18-Clean (Phthalate-Free)
 - SolEx C18-525
 - SolEx Unbonded Silica
 - Carbon Based SPE Cartridges
 - SolEx Activated Carbon
 - SolEx Graphitized Carbon
 - Polymeric New SolEx Phases
 - SolEx HRPHS
 - SolEx SAX
 - SolEx SCX
 - SolEx WAX
 - SolEx WCX



ThermoFisher
SCIENTIFIC

Rocket Evaporator



Rocket Evaporator



ASE Pucks



Flip Flop Vials

- Fully Automated Workflow: Fills gap not addressed by ASE or AutoTrace
- Reduces Laboratory Error: Pucks allow direct sample transfer from the ASE
- Expedites Sample Processing: Flip Flops evaporator directly into autosampler vials

ThermoFisher
SCIENTIFIC

Total Workflow Solution: Our Key Advantage

Key Applications:

- Fracking: ASE and LC/MS
- Cr(VI) in soil: ASE and IC
- Pesticides in beverages: AutoTrace and GC/MS

Key Regulations:

- China GB5749-2006 and GB3878-2002: SVOC in drinking water and wastewater
- EPA Method 3545A: POPs in solid matrices
- EU Water Framework Directives

Productivity:

- Thermo Scientific as the single sources of complete analytical workflow
- Reduced solvent use by as much as 90%
- Reduced extraction time to 12 min/sample



ASE 150 and ASE 350



AutoTrace 280



Rocket Evaporator



GC and GC/MS



LC/MS



IC

ThermoFisher
SCIENTIFIC

Automation of Sample Preparation in Ion Chromatography

ThermoFisher
SCIENTIFIC

pH-Conductivity Accessory

- A small unit mounts inside an AS-AP autosampler
- Conductivity detector and pH flow cells in series
 - Measurements done prior (!) to injection
 - Conditionals in CM
- USB connection for power and control/communication
- Fluidics connections
- Conductivity mode only- (no pH) just leave plug in place



ThermoFisher
SCIENTIFIC


pH-Conductivity Accessory

- Easy to install
 - USB Connection for Power and Control
- Leave plug for conductivity only
 - Install pH electrode, if pH measurement is desired



ThermoFisher
SCIENTIFIC

Conductivity based Method Selection



Column: Thermo Scientific™ Dionex™ IonPac™ CS12A,
3 × 150 mm

Eluent Source: EGC-MSA

Eluent: 20 mM MSA

Column Temp.: Ambient

Inj. Volume: A: 2.5 µL
B: 25 µL

Flow Rate: 0.5 mL/min

Detection: Suppressed conductivity

Sample: 6 cation standard in 100 mM HCl


1. Lithium	0.5	mg/mL
2. Sodium	2	
3. Ammonium	2.5	
4. Potassium	5	
5. Magnesium	2.5	
6. Calcium	5	

ThermoFisher
SCIENTIFIC

Inline Filtration with Thermo Scientific Dionex AS-DV


- Thermo Scientific™ Dionex™ AS-DV positive displacement autosampler

- Each sample uses a separate filter, NO carryover
- Filter caps (20 µm pore size)
- Very low dead volume
- Particulates tend to settle to the bottom, filter moves from top down
- Loop and concentrator loading
- No sample waste
- Low cost autosampler
- No sample pump needed



ThermoFisher
SCIENTIFIC

Dionex AS-DV: How Sample Injection Works





ThermoFisher
SCIENTIFIC

Inline Filtration with Cartridges and Valves


High pressure In-Line filters

- **Multiple configurations available**
 - Filter in sample line before injection valve
 - Inline on high pressure side
 - Dual filters and 10-port valve with backflush pump
 - High capacity, high volume filter stack of 35 μm , 7 μm and 0.45 μm filters
 - Low capacity, low volume filter frit, 0.5 μm
- **Advantages**
 - Entire sample is filtered- no sample loss
 - Small sample volume loaded onto filter maximizes capacity
 - Filter backflush minimizes carryover and extends filter lifetime



ThermoFisher
SCIENTIFIC


Inline Filtration: Easy to implement!



System Backpressure Indicates When to Exchange Filter


ThermoFisher
SCIENTIFIC

Inline Filtration: Low Pressure or High Pressure



ThermoFisher
SCIENTIFIC

Auxiliary Valves, ICS-5000 and ICS-2100



ThermoFisher
SCIENTIFIC

Thermo Scientific Dionex InGuard Cartridges

- Cartridges are high capacity

Cartridge Type	Cartridge Capacity
Ag	5 to 5.5 mEq
H	5 to 5.5 mEq
Na	5 to 5.5 mEq
HRP	2 g resin
Na/HRP	50% Na / 50% HRP



- Inline approach- sample volume per cycle equals injection volume
 - Small amount loaded onto cartridge increases effective capacity to many samples per cartridge
 - Reduced sample processing time, labor and cost

High Capacity and Small Sample Volume: Many Samples / Cartridge


ThermoFisher
SCIENTIFIC

Description of Matrix-Elimination using Thermo Scientific Dionex InGuard Products

- An autosampler or pump is used to pass the sample through the cartridge to remove matrix components.
- The cartridges have very low backpressure but the chemistries usually function optimally at flow rates between 0.5 and 2 mL/min.
- This Sample Preparation is suitable for both standard bore and microbore IC systems
- The number of samples that can be treated with an Dionex InGuard product depends on the matrix composition, concentration of compounds to be removed and the sample injection volume.

ThermoFisher
SCIENTIFIC

What is the General Setup?




Autosampler: Thermo Scientific™ Dionex™ AS-DV, AS-AP

ThermoFisher
SCIENTIFIC

Thermo Scientific Dionex InGuard Cartridges

Application example: Anion analysis in 1.6% NaCl, treat with Dionex InGuard Ag followed by Dionex InGuard Na



Column: Thermo Scientific™ Dionex™ IonPac™ AG15 and AS15, 4 mm
 Eluent: 23 mM KOH, EGC
 Flow rate: 1 mL/min
 Sample Prep: Dionex InGuard Ag and Na cartridges
 Injection vol: 100 µL
 Detection: Suppressed conductivity, Thermo Scientific™ Dionex™ ASRS 300

--- 1.6% NaCl blank without InGuard
 — 1.6% NaCl with Dionex InGuard Ag and Na
 — Water blank with Dionex InGuard
 — Std of 2ppm nitrite, sulfate, nitrate in 1.6% NaCl using Dionex™ InGuard™


Peaks: mg/L (in standard)

1	Chloride	
2	Nitrite	2
3	Carbonate	
4	Nitrate	2
5	Sulfate	2

Dionex InGuard easily removes interfering chloride, automatically

ThermoFisher
SCIENTIFIC


Anions in 50% NaOH



Column: Thermo Scientific™ Dionex™ IonPac™ AS19 w/ guard
Eluent: 15 mmol/L KOH
Transfer Eluent: Water
Flow: 0.25 mL/min
Transfer Flow: 0.5 mL/min
Sample Pretreatment: InGuard-H
Detection: Suppressed conductivity
Background Signal: < 1.5 µS/cm
Suppressor: Thermo Scientific™ Dionex™ Anion ASRS 300, recycle mode
Suppressor Current: 10 mA
Temperature: 30°C
Injection Volume: 50 µL
Sample Transfer: 1 mL
Loop:
Sample Preparation: The samples "50% NaOH" were diluted 1 ad 50 with ultrapure water and immediately analysed.


ThermoFisher
SCIENTIFIC

Matrix-Elimination Using Analyte Trapping



ThermoFisher
SCIENTIFIC

Determination of Anions in 2-Propanol




23488

25

ThermoFisher
SCIENTIFIC

Trace Anion Analysis of Electronic Grade Acetone




13508

26


ThermoFisher
SCIENTIFIC

Anions from Water Insoluble Drug




Anions Trapped and Drug Substance Removed

Instrument: Thermo Scientific™ Dionex™ ICS-3000 and AS
 Concentrator: UTAC-ULP1 (5 x 23 mm)
 Column: Thermo Scientific™ Dionex™ IonPac™ AG15 (mm)/AS15 (2 mm)
 RFIC-EG: KOH Gradient
 Temperature: 30°C
 Flow Rate: 0.4 mL/min
 Inj. Vol.: 100 μL
 Detection: Suppressed conductivity
 Sample: Drug formulation in Methanol
 Procedure:
 1. Sample passed through UTAC
 2. Anions collected; drug passed to waste
 3. UTAC flushed with water from AS to remove MeOH
 4. Anions injected onto column




ThermoFisher
SCIENTIFIC

Determination of Trace Anions in Organic Solvent



ThermoFisher
SCIENTIFIC

Determination of Trace Anions in Organic Solvent




- Collector: IonSwift™ MAC-200 (0.75×80 mm)
- Column: IonPac™ AG18-Fast capillary Guard(0.4×35mm)+IonPac™ AG18-Fast capillary(0.4×150mm)
- Switching-time: 0-1 min
- Eluent (KOH): 0-6min 5mM; 6-25min 5-25mM; 25-30min 25-45mM; 30-40min 45mM
- Flow-rate : 10µL/min
- Sample volume: 50 µL


	F	Cl	NO2	NO3	SO4	PO4
Methanol	4.7748	3.6944	13.4431	17.972	4.7748	N.D.
Acetonitrile	3.5126	2.5040	14.3851	17.7439	3.5126	N.D.
Isopropanol	3.3276	2.7701	13.9012	17.0331	3.3276	1.4484
Acetone	3.1213	2.3065	14.0539	17.1718	3.1213	N.D.

ThermoFisher
SCIENTIFIC

Anions in Hydrogen Peroxide



Concentrator Column: TAC-LP1
 Trap Column: ATC-1
 Column: Thermo Scientific™ Dionex™ IonPac™ AG15, AS15
 Eluent: KOH (EG40)
 Temperature: 35° C
 Flow Rate: 2 mL/min
 Injection Volume: 400 µL
 Detection: Suppressed conductivity




Peaks:


1. Fluoride	9.28 µL (ppb)
2. Acetate	—
3. Formate	—
4. Chloride	1.66
5. Carbonate	—
6. Sulfate	3.63
7. Oxalate	—
8. Nitrate	3.22
9. Phosphate	3.64

ThermoFisher
SCIENTIFIC

Cations in "Electronic Grade" Hydrogen Peroxide



Column: Thermo Scientific™ Dionex™ IonPac™ CS10 with Guard (4mm)
 Concentrator: Thermo™ Scientific™ Dionex™ IonPac™ CG12A (4mm)
 Eluent: 40 mmol/L MSA
 Flow Rate: 1 mL/min
 Injection volume: 1000 µL
 Detection: Suppressed conductivity
 Suppressor: Thermo Scientific™ Dionex™ CSRS-Ultra, ext. Water




ThermoFisher
SCIENTIFIC

Trace Anions in High Purity 24.5% HF

Pretreatment Column: Thermo Scientific™ Dionex™ IonPac™ ICE-AS6
 ICE Eluent: Deionized water
 ICE Flow Rate: 0.55 mL/min
 Sample Volume: 750 µL

Column: Thermo Scientific™ Dionex™ IonPac™ AG9-HC, AS9-HC, 2 mm
 Concentrator: Dionex™ IonPac™ AG9-HC, 4 mm

IC Eluent: 8.0 mM Sodium carbonate
 1.5 mM Sodium hydroxide
 IC Flow Rate: 0.25 mL/min
 Detection: Suppressed conductivity, Thermo Scientific™ Dionex™ ASRS™ ULTRA, AutoSuppression™, external water mode



Peaks:

Peak Number	Anion	Concentration (µg/L (ppb))
1	Fluoride	—
2	Chloride	7.9
3	Carbonate	—
4	Nitrate	0.89
5	Unidentified	—
6	Sulfate	10.1
7	Phosphate	2.4

14242


32

ThermoFisher
SCIENTIFIC

Trace Anions in High Purity 85% Phosphoric Acid

Ion Exclusion (ICE)
 Pretreatment Column: Thermo Scientific™ Dionex™ IonPac™ ICE-AS6
 ICE Eluent: Deionized water
 ICE Flow Rate: 0.50 mL/min
 Sample Volume: 200 µL

Column: Thermo Scientific™ Dionex™ IonPac™ AG11-HC, AS11-HC, 2 mm
 Concentrator: Dionex™ IonPac™ AG11-HC, 4 mm
 IC Eluent: 20 mM Sodium hydroxide
 IC Flow Rate: 0.38 mL/min
 Detection: Suppressed conductivity, Thermo Scientific™ Dionex™ ASRS™, AutoSuppression™, external water mode




Peaks:			µg/L (ppb)
1. Chloride	36		—
2. Unidentified	—		—
3. Sulfate	750		—
4. Nitrate	15		—
5. Phosphate	—		—

13918

33


ThermoFisher
SCIENTIFIC

Determination of Trace Anions in Weak Acids



ThermoFisher
SCIENTIFIC

Trace Anions in Formic acid – Capillary IC




1st Dimension
 Column: Thermo Scientific™ Dionex™ IonPac™ ICE-AS1(250×9 mm) ion-exclusion column
 Sample volume: 50 µL

2nd Dimension
 Concentr. Col.: Thermo Scientific™ Dionex™ IonSwift™ MAC-200 (0.75×80 mm)
 Column: Thermo Scientific™ Dionex™ IonPac™ AG18-Fast w/ guard (0.4×150mm)
 Eluent: KOH-Gradient (RFIC)
 Flow rate: 10 µL/min
 Detection: Suppressed Conductivity
 Suppressor: Thermo Scientific™ Dionex™ ACES (Recycle Mode)

ThermoFisher
 SCIENTIFIC

Trace Anions in Acetic Acid – Capillary IC



1st Dimension
 Column: Thermo Scientific™ Dionex™ IonPac™ ICE-AS1(250×9 mm) ion-exclusion column
 Sample volume: 50 µL

2nd Dimension
 Concentr. Col.: Thermo Scientific™ Dionex™ IonSwift™ MAC-200 (0.75×80 mm)
 Column: Thermo Scientific™ Dionex™ IonPac™ AG18-Fast w/ guard (0.4×150mm)
 Eluent: KOH-Gradient (RFIC)
 Flow rate: 10 µL/min
 Detection: Suppressed Conductivity
 Suppressor: Thermo Scientific™ Dionex™ ACES (Recycle Mode)

ThermoFisher
 SCIENTIFIC

Thank you!



ThermoFisher
SCIENTIFIC

The world leader in serving science

ThermoFisher
SCIENTIFIC